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OB-J FA

OB-J

	<b>MANUALE D'USO</b>	<b>3 - 13</b>
	<b>USER MANUAL</b>	<b>14 - 23</b>
	<b>BEDIENUNGSANLEITUNG</b>	<b>24 - 34</b>
	<b>MANUAL D'UTILISATION</b>	<b>35 - 45</b>

CE 0123

PRODOTTO DA:

## INDICE

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
## INDEX

INDEX	14
SYMBOLS AND IMPORTANT SAFEGUARDS	15
WARNINGS, PRECAUTIONS AND INFORMATION	16
<b><i>OB-J COLLECTION JAR</i></b>	
MODELS	17
GENERAL DESCRIPTION	17
OB-J FA COLLECTION JAR	17
PROTECTION FILTER	17
OB-J COLLECTION JAR	18
FEATURES OF THE COLLECTION JAR	18
DEVICE INSTALLATION ON MEDICAL SUCTION UNITS	19
DEVICE INSTALLATION WITH CENTRALIZED VACUUM SYSTEMS	19
CONTROLS AND CHECKS	19
USE PROCEDURES	20
MAINTENANCE AND REUSE	20
ACCESSORIES, SPARE PARTS	22
SERVICE	22
TECHNICAL DATA AND CONFORMITY TO INTERNATIONAL LAW	23
GUARANTEE	23
DECLARATION OF CONFORMITY	46

## SYMBOLS

### S.1 Symbols used on the device and recalled in the user manual

	Specific warnings regarding the device, which must always be taken into consideration.
	CE mark perform to the MDD 93/42/EEC for Medical devices with class above I
	Manufacturer
	Production date
	The materials composing the device can be recycled following specific procedures provided by national laws and local regulations.

### S.2 Symbols used in the user manual to recall the reader's attention

	Important safety information to avoid the risk of injury to the user or patient and / or damage to the device
	Information to which you should pay special attention
	Note or information for correct use of device

## CAUTION NOTICE: READ CAREFULLY!

- Please read this User Instruction manual carefully before using the OB-J FA/LINER collection jar.
- The OB-J FA/LINER collection jar is an independent medical device. Federal law restricts this device to sale, distribution and use by/or upon order of a physician, emergency medical technician, or other medical practitioner.
- Due to the nature of the jar use, all precautions must be taken to ensure the jar reliability (cleaning, storage, environmental conditions, etc.) regular operational and safety checks should be made.
- Act as indicated in chapter 4 of this manual after each use.
- If connected to the main device all suction parts need to be in compliance with ISO10079-1 (if active devices) and ISO10079-3 norms (if connected to centralized vacuum sources). Please refer to the above mentioned norms also with relation to pneumatic connection, vacuum regulators and protection filters. An eventual wrong connection of the device might cause operator's and suction circuit contamination.
- Do not use different brackets from those indicated by the manufacturer for jar fixing. Defective and/or uncorrect bracketing might damage both suction system and people using it.
- Do not use the collection jar without protection filter.
- In case of sterilization adopt all required precautions described in the proper chapter of this manual.
- It is strongly recommended that the collection jar is replaced after 30 autoclave cycles.
- There are no serviceable components in the collect container, always refer to trained service people only.
- Only use official and authorized spare parts supplied by the manufacturer (Oscar Boscarol srl company) or their representatives. Using original spare parts will grant the unit reliability and will ensure warranty validity.

### WARNING ON REUSE OF DISPOSABLE PARTS

- Reuse of disposable materials and parts may cause loss of device functionality and represent direct and indirect sources of contamination for users and patients alike.
- Sterilization and/or cleaning of consumable parts (antibacterial filters, disposable bags, suction tubes, Jankauer, etc...) can cause structural degradation of material and consequent dangers arising from loss of mechanical integrity.

## IMPORTANT INFORMATION

### Preventive maintenance and safety inspection:

The device should be inspected or replaced at least once every 24 months. Refer to the periodical checks section and always rely on the manufacturer or the authorized service centers

### Contamination of the device:

Shipment of contaminated devices is strictly forbidden. Each devices received under such conditions will be refused and health authorities will be properly informed, due to contamination risks.

## OB-J COLLECTION JAR

### Available models:


- BSU500      OB-J FA 1000 ml collection jar (without protection filter)  
 BSU506      OB-J collection jar (without disposable bag SERRES)

### General description

OB-J collection jar has been designed to be used in the medical field in order to collect and reserve organic fluid aspired. The device, connected to a vacuum source (suction units or centralized vacuum systems) is used to release patient's oropharyngeal tract so to facilitate both spontaneous and forced breathing. The device has been designed and manufactured in compliance with ISO10079-1 norm, particularly for HIGH VACUUM-HIGH FLOW applications. Easy assembling and the absence of metallic parts makes of this a resistant, lightweight and handy device.

### OB-J FA collection jar

OB-J FA model includes one jar, one lid, one over-flow valve and one 90° plastic joint. The lid of the canister includes a protection filter which is directly integrated. The OB-J FA jar can be sterilized by standard steam autoclave machines at a maximum temperature of 121 °C and at a pressure of 2 bar (200kPa). The jar should be replaced every 30 sterilization cycles. The jar must be used vertically to avoid the intervention of the over-flow valve. In case the valve is actioned, switch off the device and disconnect the bottle hose indicated with the VACUUM word indicated on the cover.

### Protection filter (separately supplied)

To prevent fluid overflow, a special protection filter is used between the jar and the unit. The filter is produced with PTFE hydrophobic material which prevents fluids from entering the pneumatic circuit. Working together with the shut-off valve on the jar, the filter isolates the pneumatic suction pump from gas and fluid substances. When the filter is wet, it is not possible to use the unit; for this reason, the filter should be changed immediately.

#### Warnings and precautions on the use of protection filter

**If the suction unit is used either in an emergency or on a patient where the risk of contamination is not known, the filter must be changed after each use.**

Even if the risk of contamination is minimal, it is recommended that a new filter is fitted each time the unit is used; otherwise, the filter should be changed monthly in any case. The material may experience structural degradation due to environmental and external conditions (humidity, heat, cold).

**Do not use the suction unit without protection filter and collection jar!**

### OB-J collection jar

OB-J includes the jar (1), adaptor for disposable bag SERRES (2) and the fitting "L" (3). The disposable bag is provided separately. The integrated filter in the bag performs the function of preventing backflow of aspirated liquid into the suction unit, when it is completely filled.

The OB-J jar can be sterilized by conventional steam autoclave at a maximum temperature of 121 °C and a pressure of 2 bar (200kPa). Every 30 sterilization cycles, the jar should be replaced.

### Features of the collection jar

#### Device installation on medical suction units operating from an electric source

OB-J collection jar may be connected to electric suction units whose aspiration ought not to exceed 950mbar (95kPa). The connection needs to be performed with a tube to the protection filter placed on the cover of the jar itself. The tube needs to have a 6mm internal diameter and a 11mm external one. Its maximum length does not have to exceed 180cm. Using transparent, strengthen and anti-narrowing tubes is preferable. Connection tubes must necessarily be separated from suction source for service and cleaning reasons.


Fig.1

Fig.2

Using different tubes from the ones indicated, might compromise the device correct functioning. OB-J collection jar has been designed and manufactured referring to actual specific industrial standards. Using tubes with an internal diameter inferior to 6mm does not assure the suction of solid substances that might occlude the tubes themselves. On the contrary, suction tubes with exceeding diameters do not allow proper connection to the jar cover causing possible loss of substances.


**NOTE**

The Oscar Boscarol Company is at disposal for eventual supplying of connection tubes peculiar to the vacuum (suction) and destined to the connection with the OB-J jar. Contact authorized distributors, assistance points and/or manufacturer for further information.

**Installation of the collection jar with centralized vacuum systems (nosocomial-like for example)**

OB-J collection jar may also be connected to centralized vacuum systems (such as nosocomial ones). In this case carefully follow the instructions as showed in the picture for a correct connection. Whether the disposable liner or the autoclavable bottle is used, we suggest to use the optional support bracket. The connection may occur with a proper tube with an internal diameter of 6mm. The connection tube maximum length does not have to exceed 180cm. The use of transparent, strengthen and anti-narrowing tubes is preferable.

Support bracket

**Controls and checks**


Check the correct functioning of the suction system connected to the OB-J jar, as per hospital routine procedures (daily, anyhow). Always check the correct holding of the bracket anchored to the wall.

The jar was designed and constructed to be safe and reliable in its operation. Removable parts have been sensibly reduced for an easy reuse, cleaning and sterilization. In case of suspected malfunction, it is better to isolate the device and contact the authorized service center. In case of failure of over-flow valve replacing the jar is necessary. The protection filter acts as a second safety protection in the event of failure of over-flow valve. Check pneumatic

connections paying attention to the connection to the central system. Only use tubes with a proper diameter (internal 6mm, external 11mm). The patient circuit connection (normally with a 130cm long tube at least) needs to be performed on the side of the cover stating "PATIENT". It is possible to connect a sterile junction (Finger-Tip) and catheter with proper dimensions to the suction tube.

The bottle should be replaced every 24 months or after every 30 sterilization cycles. Check the device as indicated below.

Fig.4

**NOTE**

If you malfunction is suspected, completely replace the jar! Always check the vertical position of the jar (cover up).

### Use procedures

Please act as follows in order to preserve reliability and security of the device, accordingly to its use (if connected to an independent suction source or to a vacuum centralized system):

- Ensure perfect capacity of jar cover before activating vacuum source.
- Check correct casing of protection filters in their proper place.
- Always check regular functioning of shut-off valve (free movement of floater).
- Ensure that control vacuum level can be easily checked.
- During suctioning, always verify possible intervention of shut-off valve.
- Do not connect the jar to suction sources with negative pressure superior to 95kPa (712mmHg or 950mbar)

#### WARNING

**In case of intervention of shut-off valve, immediately arrest suction source. Do not leave the jar connected to the suction device after valve intervention for a period of more than 5 minutes.**

## MAINTENANCE AND REUSE

### Occupational safety and PPE (Current legislation concerning workplace safety and health)

The decontamination process is always a delicate process, which implies a specific training, especially in the field of emergency, where the patient's medical condition and contamination level is mostly unknown. For this reason, the user must always use personal protective equipment (PPE) to protect himself and other people.

In case a proper PPE is not available, please contact the person who is in charge for the workers' safety in your company or organisation.

The organic secretions collected by the device's jar might cause serious infections to the users. For this reason, suitable PPE and disinfection products must always be used, applying what has been established by the competent authorities.


### After each use of the OB-J FA collection jar

After each use, switch off the suction unit and disconnect the connection tube.

Cleaning and disinfecting of the reusable jar:

- Always wear gloves and protective clothing.
- Carefully remove the jar from the unit.
- Remove yellow angular connector with patient tube (2). If the tube is equipped with a Yankauer probe, it must be disposed of together with the curved probe (disposable devices). The angular connector can be sterilized.
- Remove the protection filter from the jar (3).
- Remove the red cover (4). Attention to possible contamination with the liquids contents! Empty the contents of the jar in accordance with local regulations or hospital practice.
- Empty the jar and dispose of the contents, the protection filter, the disposable Finger-tip and the suction catheter in accordance with local regulations or hospital practice.
- Disassemble the jar (5-8).

During assembly and disassembly always check that jar parts are free from defects and malfunctions. The over-flow valve includes a float that slides into a plastic lead. Make sure of its correct sliding. To remove filter from cover, respectively unscrew or screw it in its housing. This operation facilitates the extraction and avoids possible damage the cover.

### Jar decontamination

If sterilization is necessary, only use steam autoclaves at a maximum temperature of 121°C for 15 minutes. Do not use pressure exceeding 2bars. The canister of the jar should be placed in the autoclave upside down (bottom up). After the sterilization cycle cool down to environmental temperature, check its integrity and reassemble the bottle following the inverse operations used to dismantle it.

#### CAUTION:

- Do not load weight to the parts during sterilisation cycle.
- Observe maximum temperature limits, pressure and duration during the autoclave cycle.
- Never exceed the value of 60°C for washing or cleaning operations (with the exception of sterilization in a steam autoclave).
- Cleaning, disinfection and sterilisation should only be carried out by trained personnel.
- The collection jar can be used up to 30 sterilisation cycles, after that jar, tubes and plastic connectors must be replaced.
- After jar reassembling, check that the lid is properly fitted so to avoid substance lacks.

### After each use of the OB-J collection jar




After each use, remove the collection jar (with the liner inside) from the unit and put it on a flat, level surface. Always wear gloves and protective clothing when replacing the collection jar liner. Disconnect the patient tube from the white connector on the disposable bag (2) and close the connection <PATIENT> with the plug provided on the lid (3). Remove the sealed liner from the jar (4). Dispose of the liner in accordance with local regulations or hospital practices. Disconnect the silicone tube from “L” connector (5). Remove the adaptor for disposable bag (6) and unscrew the “L” connector.

Soak the parts in cold water and rinse well. Wash all parts of the jar with hot water (max. temperature 60°C). After the pre-washing with water, it is possible to use substances, disinfectant and/or bactericides, avoiding the use of substances, which can irreparably damage the jar (for example Betadine). After disinfection, rinse and dry all parts. If it is necessary to sterilize the jar, proceed as described on the previous page “Jar decontamination”.

### Reassembly of the jar

Extract a new disposable bag from the package, stretch it (picture 1) and insert it into the jar (picture 2). Connect the complete jar to the suction unit. Activate the suction unit. Close the connector <PATIENT> with a finger and, at the same time, press the bag lightly from the center of the lid (picture 3). Make sure the bag is completely swollen. Connect the patient tube (Jankauer) to the connector <PATIENT>.

**The disposable bag must be replaced after each use!**

### Disinfection and decontamination of jars

To correctly disinfect and decontaminate the suction unit, we recommend using specific, approved products. These disinfectants must be free of alcoholic and/or abrasive substances. Oscar Boscarol srl can provide specific materials for disinfection of medical equipment, including our suction units. These disinfectants, available in different formats (wipes, spray, liquids), have been laboratory tested and guaranteed to deactivate viruses, bacteria and microorganisms. When used periodically, they destroy and prevent the formation of dangerous biofilms (superficial layers that easily host bacteria, moulds, viruses and microorganisms). Our disinfectants do not contain alcohol, chlorine, phenols, aldehydes and halogens.

**NOTE**

For more detailed information, contact us at [info@boscarol.it](mailto:info@boscarol.it) visit our website [www.boscarol.it](http://www.boscarol.it).

## ACCESSORIES AND SPARE PARTS

Index code	Description
<b>Parti di consumo</b>	
BSU730	Protection filter for OB-J FA jar – 5 pcs
BSU732	Protection filter for OB-J FA jar – 15 pcs
BSU734	Protection filter for OB-J FA jar – 40 pcs
BSU705	Disposable bag SERRES – 6 pcs
BSU706	Disposable bag SERRES – 12 pcs
BSU707	Disposable bag SERRES – 36 pcs
<b>Parti di ricambio</b>	
SPS6000	Bottle OB-J FA without lid and without protection filter
SPS6002	Shut-off valve – 3 pcs
SPS6004	90° plastic joint for OB-J FA jar – 3 pcs
SPS6006	Lid for OB-J FA complete with shut-off valve and 90° plastic joint
SPS5090	Bottle OB-J complete, without disposable bag
SPS5092	“L” fitting - 3 pcs

**NOTE**

The parts listed above may be changed without prior notice from the manufacturer for technical improvement reasons. Contact the manufacturer for more information.

## SERVICE

The device is do not need any repair. In case of plastic parts faults they must be replaced.

## TECHNICAL DATA AND CONFORMITY TO INTERNATIONAL LAW

### Classification according to the MDD 93/42/EEC

Collection jar is in compliance with ISO10079-1 for the connection to the active suction unit and ISO10079-3 norms if connected to centralized vacuum sources.

MDD Classification: **Ila**

### Dimensions

Maximum diameter of OB-J FA:	ø 106 mm
Minimum diameter of OB-J FA:	ø 88 mm
Maximum height of OB-J FA:	230 mm
Weight of OB-J FA (without protection filter):	168gr. ± 5gr
Maximum width of OB-J (without disposable bag):	125 mm
Minimum diameter of OB-J:	ø 88 mm
Maximum height of OB-J:	210 mm
Weight OB-J (without disposable bag):	176gr. ± 5gr
Tolerance on all values:	±5%

### Technical data

Maximum vacuum applicable:	950mbar (95kPa), 712mmHg +/-10% on the whole range for max.5 minutes
Jar gradual change range interspaced:	50ml ± 5%
Max. capacity:	70 L/M
Type of material used for fabrication:	RG835MO (polypropylene for medical and pharmaceutical appliance).
Sterilisation of the collection jar:	121°C for 15 minutes at per hospital procedures (at a maximum pressure of 2bar)
Mechanical securities	Shut-off valve made from polypropylene
Connection tubes on patient side:	6mmx9mm (max. length 180cm)
Connection tubes on vacuum side:	6mmx9mm (max. length 200cm)

### Environmental and transport conditions

Functioning temperature	0÷50°C
Storage temperatures:	-20÷60°C
Relative humidity for running or storage:	20÷80%
Range of atmospheric pressure:	70÷130kPa (700÷1300mbar; 525÷975mmHg)

## GUARANTEE

The Oscar Boscarol Company guarantees the OB-J collection jar for a period of 12 months from the date of purchase by original user. Oscar Boscarol srl guarantees the absence of defects in material and workmanship on the suction unit. Unauthorized technical interventions, modifications, mistaken connection to the suction circuit and everything which does not relate to the device intended use is excluded from this warranty (including **wear and tear of collection jar, discoloration, and any other aesthetic irregularities which do not affect the unit functionality**). All shipping costs are borne by the customer. With reference to the conditions described in this warranty, Oscar Boscarol srl. cannot be held liable for incidental or consequential damages due to modifications, repairs, unauthorized technical intervention or whether one of its part had been accidentally damaged by misuse or abuse of any sort. No other expressed or limited warranties are provided on the device concerning commercial or suitability aspects but those described in this user manual. For any arising legal matters the Court of Bolzano (Italy) is the one in charge.

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