



BSU216

OB1000 FA MEDICAL SUCTION UNIT + EU SUP.

- Medical Suction Unit OB 1000 powerful, easy and compact
- Max. suction flow 25 liter/minute
- Max suction rate 850 mbar (85 kPa)
- Autonomy on continuous functioning >45 Minutes.
- Autoclavable 1000 ml jar with safety valve and protection filter integrated on the lid
- Cable for the connection to the vehicle lighter (12Vcc)
- Provided with a 2-poles male connector main power supply, (power supply voltage rate between 100 e 240 Vac 50/60 Hz)
- Dimensions: about 34x11x24 cm
- Weight: about 3,3 kg
- Manufactured in compliance with all current norms

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.